MiniLift200 SystemRoMedic





A unique solution for active sit-to-stand - with electrical base widening

MiniLift200 is a mobile sit-to-stand lift which is used, in combination with a specially adapted lifting vest, to assist users with reduced function when rising from a sitting to a standing position. MiniLift200 features a construction offering a totally unique and extremely comfortable sit-to-stand experience. The lift is also compact and convenient and very easy to manoeuvre; features which, in combination with the possibility to adjust the settings for an optimal fit for every user, makes it very useful in most locations and situations.





MiniLift20

Extremely comfortable, active sit-to-stand with natural movement pattern

MiniLift200 has been developed to, as gently and comfortably as possible, assist the user when rising from a sitting to a standing position. The handles and the side supports offer multiple hand placement options for the user. The stable anti-slip foot plate is very low and the lower leg supports, as well as the lifting arm, can be adjusted to fit every user optimally. When MiniLift200 is used in combination with the appropriate lifting accessories, the user is supported under the feet, for the front of the lower legs and behind the back, which provides for a safe and secure, yet active sit-to-stand procedure. The construction safely moves the user forwards and upwards in a natural movement pattern and, at the same time, the leg muscles and balance are exercised. MiniLift200 is manufactured in steel for stability and strength and is suitable also for tall and heavy users weighing up to 200 kg/440 lbs.

Active sit-to-stand with natural movement pattern.



PRODUCT SPECIFIKATION

MiniLift200 __Mobile sit-to-stand lift with electrical base widening _ Art. No. 401100334

TECHNICAL INFORMATION

Safe working load (SWL)	200 kg / 440 lbs
User length, min-max	1400-2000 mm / 55.1-78.7"
Material	Powder-coated steel
Height	1000 mm / 39.4"
Length	980 mm / 38.6"
Width of base, outer measurement	590-1030 mm / 23.2-40.6"
Width of base, inner measurement	450-875 mm / 17.7-34.4"
Weight	49 kg / 108 lbs
Base widening	Electrical
Height of footplate, for user	95 mm / 3.7"
Castor size, diameter	100 mm / 3.9"
Lifting speed, without load	35 mm / 1.4" per sec
Emergency lowering	Manual and electrical



MiniLift200 meets the requirements of the Medical Devices Directive for Class 1 products 93/42EEC.





New control box with diagnostics, e.g. lift counter and overload

Electrical emergency stop button readily available on the control

ACCESSORIES







45690004-008

Art. No.	Product name	Description
45500003 45500004 45500006 45500007 45500008 45500009	ThoraxSling ThoraxSling ThoraxSling ThoraxSling ThoraxSling ThoraxSling ThoraxSling	Polyester, X-Small Polyester, Small Polyester, Medium Polyester, Large Polyester, X-Large Polyester, XX-Large
45550004 45550006 45550007 45550008 45550009	ThoraxSling ThoraxSling ThoraxSling ThoraxSling ThoraxSling	Wipeable, Small Wipeable, Medium Wipeable, Large Wipeable, X-Large Wipeable, XX-Large
45590004 45590006 45590007 45590008	ThoraxSling ThoraxSling ThoraxSling ThoraxSling	Disposable, non-woven, Small, 5 pcs Disposable, non-woven, Medium, 5 pcs Disposable, non-woven, Large, 5 pcs Disposable, non-woven, X-Large, 5 pcs
45600004 45600006 45600007 45600008	ThoraxSling ThoraxSling ThoraxSling ThoraxSling	With seat support, polyester, Small With seat support, polyester, Medium With seat support, polyester, Large With seat support, polyester, X-Large

ThoraxSling With seat support, disposable, non-woven, Small, 5 pcs With seat support, disposable, non-woven, Medium, 5 pcs With seat support, disposable, non-woven, Large, 5 pcs ThoraxSling ThoraxSlino

With seat support, disposable, non-woven, X-Large, 5 pcs For MiniLift200 as from serial number 200001277 Soft, wipeable protection for the bar of MiniLift200 6 buttons, with service, battery status and overload indicators External charger for all models and variants of MiniLift

ThoraxSling

ThoraxSling is a specially adapted and very comfortable sit-to-stand vest that is used in combination with Handicare's SystemRoMedic™ MiniLift sit-to-stand lifts. It is easy for the caregiver to apply, while giving the user an added sense of security during rising and standing training. ThoraxSling also facilitates toileting. ThoraxSling is available in two models; one which provides support behind the back (ThoraxSling) and one that also gives support under the seat (ThoraxSling, with seat support), and in several materials and sizes.

Other models and variants of MiniLift mobile sit-to-stand lifts

MiniLift is available in several models and variants with different lifting capacity, with electrical or manual base-width adjustment and with standard height or low legs for use with lower beds.

Art. No.	Product name	Description
60300010	MiniLift160EM	Mobile sit-to-stand lift with manual base widening, SWL 160 kg/350 lbs
60300011	MiniLift160EM	Mobile sit-to-stand lift with manual base widening, low legs, SWL 160 kg/350 lbs
60300012	MiniLift160EE	Mobile sit-to-stand lift with electrical base widening, SWL 160 kg/350 lbs
60300013	MiniLift160EE	Mobile sit-to-stand lift with electrical base widening, low legs, SWL 160 kg/350 lbs
400641434	MiniLift125	Mobile sit-to-stand lift with electrical base widening, SWL 125 kg/275 lbs

Please visit www.handicare.com/us for more information



Toll Free: 866 276 LIFT/5438 www.handicare.com/us

customerservice.usa@handicare.com

45690004

45690006

45690007

45690008

70200033

70200012

ThoraxSling

SlingBarWrap MiniLift

Hand control HB33-6 External charger Linak

CalfStrap