

The Jazzy[®] Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.

> Jazzy® Elite HD shown in Jazzy Red with high-back seat



SPECIFICATIONS

Weight Capacity	450 lbs.*
Maximum Speed	Up to 4 mph
Ground Clearance	3"
Turning Radius	24"
Overall Length	35.5"
Overall Width	24.25"
Drive Wheels	14″ solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, front-wheel drive
Braking System	Intelligent braking (electronic, regenerative disc brakes)
Standard Electronics	60 amp, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 9.8 miles
Battery Requirements	(2) 12 volt, deep cycle, NF-22
Battery Weight	37 lbs. each
Base Weight	98.5 lbs.
Standard Seat Weight	51.5 lbs. (high-back)
Standard Seat	20" x 18-20"
Max. Seat Dimensions	24" x 24"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors



OPTIONAL SEATING

- Synergy[®] Seat
- Solid seat pan

MODELS	
Jazzy Elite HD-2HD-SS	Code: K0824
Jazzy Elite HD-2HD-C	Code: K0825

FEATURES

- · Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 70" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag
- Heel loops

- Lap belt (50", 60", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Residual limb support
- Military patches
- Swing-away joystick
- XLR USB charger

RECOMMENDED LIFTS

- Backpacker[®] Plus
- Backpacker[®] AVP
- Outlander Full Platform
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com



(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or – 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.